



A FUSION OF COMFORT AND INNOVATION

A Fusion of Comfort and Innovation

The Hybrid is a unique fusion of technologies, designed to be **a universal interface for patients, clinicians and providers.** This one product serves as a nasal interface, and an oral/nasal interface.

The **Hybrid is a dual airway interface** for use with positive pressure devices for the treatment of respiratory insufficiencies and obstructive sleep apnea.

The unique design of the Hybrid promotes patient comfort by **eliminating pressure points** on the forehead and across the bridge of the nose.

The Hybrid complete system **contains all of the components required to properly size the patient**, eliminating the hassle of opening multiple packages to get the right fit.



Hybrid Universal Interface

Performance features

- Built-in chin flap gently supports the chin, eliminating the need for a chin strap
- Premium headgear provides five points on adjustability for better fit, seal and increased support
 - Patient can wear glasses while wearing the interface
- Quick connect/disconnect on headgear enables easy removal and application
 - No forehead or overhead contact

Design Features

- 100% latex free
- One package contains all of the components required to properly size the patient, eliminating the hassle and waste of opening multiple products to properly fit the patient
- Universal design simplifies choice and fitting protocols and minimizes time spent with a patient to address issues such as mouth breathing or claustrophobia
- Adjustable nasal pillow design provides a customized fit for each patient
- No nasal bridge or forehead contact eliminates nasal bridge sores and air leaks near eyes
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports
- Patient Instructional CD and Carrying case included with each complete system
- To view video, please visit : www.hybridmask.com

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
HYBRID MASK (HYB500)	Combination oral/nasal mask, used with continuous positive airway pressure device. each	A7027	1 every 3 months	\$179.35-\$186.52
HYBRID ORAL CUSHION (HYB511, HYB513, HYB515)	Oral cushion for combination oral/nasal mask, replacement only. each	A7028	2 pairs per month	\$49.54 -\$49.54
HYBRID NASAL PILLOW (HYB521, HYB523, HYB525)	Nasal pillows for combination oral/nasal, replacement only. pair	A7029	2 pairs per month	\$20.24 -\$20.24
HYBRID HEADGEAR (HYB530)	Headgear	A7035	One per 6 months	\$33.75 - \$33.79

Accessories & Replacement Parts

RESPCARE

PART NUMBER

PRODUCT DESCRIPTION

HYB500 Hybrid Interface Complete System
 HYB515 Hybrid Interface Oral Cushion Large
 HYB513 Hybrid Interface Oral Cushion Medium
 HYB511 Hybrid Interface Oral Cushion Small

PART NUMBER

PRODUCT DESCRIPTION

HYB530 Hybrid Interface Headgear
 HYB525 Hybrid Interface Nasal Pillows Large
 HYB523 Hybrid Interface Nasal Pillows Medium
 HYB521 Hybrid Interface Nasal Pillows Small

09030 REV A

6601 Lyons Rd. Suites B1-B4, Coconut Creek, FL 33073

Order : 888-280-6578