

Humidifier/device compatibility list

This supplement specifies the devices compatible with ResMed's humidifiers. For full details regarding the correct use of these devices, please refer to the relevant User Guides.

Device \ Humidifier	Passive Humidifiers		Heated Humidifiers			
	Passover	HumidAire 2iC™	HumidAire 2i™	HumidAire 3i™ H3i Plus™	H4i™	H5i™
AutoSet CS™ 2		✓	✓			
AutoSet Spirit™	✓	✓	✓			
S7™ Series	✓	✓	✓			
S8 AutoScore™ S8 AutoSet Vantage™ S8 AutoSet Spirit™				✓	✓	
S8 Escape™ S8 Lightweight™	✓			✓	✓	
S8 Escape™ II				✓	✓	
S8 Elite™ II S8 AutoSet™ II S8 AutoSet Spirit™ II				✓	✓	
S9™ Series						✓
VPAP™ III Series	✓	✓	✓			
VPAP™ III ST-A with QuickNav		✓	✓			
VPAP™ Adapt SV		✓	✓			
VPAP™ Auto				✓	✓	
VPAP™ Auto 25 VPAP™ ST VPAP™ IV					✓	
VPAP Malibu™			✓			

Note: Not all devices are available in all regions.