

## Mirage Swift™ Nasal Pillows System

### COMPONENTS CARD

#### Fully Assembled Mirage Swift Nasal Pillows System

Nasal Pillows System **60505**  
(Includes a small, medium and large pillow sleeve)

Frame System with pillow sleeve  
(excluding headgear)

- 60532** (small)
- 60533** (medium)
- 60534** (large)

Swift Bag  
**60531**

Frame  
(with clip and cap)  
**60527**

Elbow

Short tube assembly  
(elbow, short tube  
and swivel)  
**60528** (1 pack)  
**60529** (10 pack)

Tube retainer  
**60530**  
(3 pack)

Swivel

Clip  
**60535**

Distal pressure port  
**61121**  
(not part of the system.  
Can be ordered separately)

Headgear buckle  
**60523**  
(4 pack)

Headgear  
**60526**

Seal ring  
**60524** (10 pack)

Pillow sleeve  
**60520** (small)  
**60521** (medium)  
**60522** (large)

Cap  
**60525**

Specifications may change without notice.  
Mirage and Mirage Swift are trademarks of ResMed Ltd.  
©2004 ResMed Ltd. Reference number: 1010220/2 04 12

